

MIKRO 200/200 R Microlitre Centrifuges classic/cooled



# MIKRO 200/200 R

The MIKRO 200 and 200 R rank among the fastest centrifuges in their class – with a maximum speed of 15,000 RPM and an RCF of 21,382. Not only are they fast, the MIKRO 200 R also refrigerates quickly.

Even thermosensitive samples can be gently centrifuged thanks to highly reliable refrigeration. Best separation results are guaranteed with the smooth-running, high-performance motor featuring a maintenance-free frequency drive.

With the autoclavable and aerosoltight accessories, even infectious materials can be processed safely for both the user and the environment.

## **USER-FRIENDLINESS**



- Powered lid locking
- Quick-entry foil keypad
  - · Easy-to-read, digital display
  - · Display of current parameter values
  - · Select key for parameter selection
  - · Impulse key for short centrifuging
  - · Open lid key
- · RCF key
- Twist knob to enter values
  - · RPM in increments of 10
  - · Running time in minutes, max. 99 min
  - Temperature in increments of 1°C (MIKRO 200R)
- 4 programmable memories

# REFRIGERATION (MIKRO 200 R)

- Infinitely variable setting from -10 °C to +40 °C
- Fast Cool function
- Stand-by mode maintains the set temperature
- Automatic refrigeration switch-off when the lid is opened

#### **SAFETY**

- Lid dropping protection
- Lid locking and holding
- Emergency lid lock release
- Motor overheating protection
- Chamber overheating protection
- Imbalance switch-off
- Standstill indication
  - with the refrigerated model MIKRO 200 R the LED in the Open key lights up
  - with the MIKRO 200 the lid opens automatically after the rotor has come to a standstill



## **OPERATION AND CONTROL**

The display panel and all controls on the MIKRO 200 and 200 R centrifuges have been designed ergonomically. Information on the display is clearly positioned and perfectly legible. During centrifugation, the actual parameter values are indicated.

Nine acceleration and braking ramps, or unbraked run-down, can be preselected. Ramp 9 always corresponds to the shortest possible run-up or run-down time.

The stored parameter combination remains in the memory even after the centrifuge has been switched off.

By entering the rotor radius r/mm, the centrifuge calculates and indicates the rotational speed RPM or the relative centrifugal force RCF.

The twist knob allows quick and easy entry or changing of the parameters.

## **KEYPAD**

Precooling (MIKRO 200 R)

RCF Change-over from RPM to RCF indication

SELECT Selection of parameters

START Starts centrifugation /
IMPULS Short cycle mode

STOP Stops centrifugation manually /

OPEN Opens the lid

## **ENTRY OF THE PARAMETERS**

P Entry of the program number,

4 programmable memories

T/°C Entry of the temperature in increments

of 1 °C within a range of – 10 °C to + 40 °C

(MIKRO 200 R)

RCF Entry in increments of 10 RPM Entry in increments of 10

r/mm Entry of the rotor radius in mm

t/min Entry of the centrifuging time in minutes and seconds (max. 99 min : 59 sec)

Entry of the acceleration ramp 1 - 9 Entry of the braking ramp 1 - 9



Keypad of the MIKRO 200 R (N Plus control panel)



# **ROTORS AND ACCESSORIES**



capacity in ml	0.2	0.4	0.5	0.8	1.5	2.0
Ø x L in mm	6 x 18 6 x 45 8 x 3		8 x 30	8 x 45	11 x 38	
Cat. No.			-			0536
lid <b>E3243</b> incl.	Į					
rotor Cat. No. 2424-B			9			
Cat. No.	2024		2023		<b>2031</b> <sup>2)</sup>	-
boring Ø x L in mm	6 x 40		8 x 40		10.2 x 19	11.2x42.6
tubes per rotor	24					
max. RCF	21,382					
radius in mm	85					
run-up in sec	20					
run-down in sec, braked	28					
temperature in °C1)	+4					

capacity in ml	0.2	0.4	0.5	0.8	1.5	2.0
Ø x L in mm	6x18	6x45	8x30 8x45		11x38	
Cat. No.	-				2078	0536
lid E3243 incl.		9				
rotor Cat. No. 2427-A						
Cat. No.	2024		2023		<b>2031</b> <sup>2)</sup>	-
boring Ø x L in mm	6×40		8×40		10.2×19	11.2 x 41.3
tubes per rotor	30					
max. RCF	20,627					
radius in mm	82					
run-up in sec	22					
run-down in sec, braked	30					
temperature in °C1)	+4					

capacity in ml	0.5
Ø x L in mm	8 x 30
Cat. No.	-
lid <b>E3243</b> incl.	
rotor Cat. No. 2430-B	
	-
Cat. No. <b>2430-B</b>	- 8.1 x 30.3
Cat. No. 2430-B	- 8.1 x 30.3 30
Cat. No. 2430-B  Cat. No.  boring Ø x L in mm	
Cat. No. 2430-B  Cat. No.  boring Ø x L in mm  tubes per rotor	30
Cat. No. 2430-B  Cat. No.  boring Ø x L in mm  tubes per rotor  max. RCF	30 20,376
Cat. No. 2430-B  Cat. No.  boring Ø x L in mm  tubes per rotor  max. RCF  radius in mm	30 20,376 81

<sup>1)</sup> Lowest attainable temperature in precooled refrigerated centrifuges at max. speed. Lower temperatures can be attained by reducing the speed. 2) For centrifugation at high speeds, we recommend to use form-fitting, phenol-resistant adapters 2031.

## **ROTORS AND ACCESSORIES**

## Angle rotor, 20-place, for cryo tubes

# ₹40°

n = 15,000 min<sup>-1</sup> max. RCF 21,130

Cat. No. 2426-B

## Angle rotor, 4-place, for PCR strips



**⋠**45°

 $n = 15,000 \text{ min}^{-1}$  max. RCF 14,338

Cat. No. 2418-A

capacity in ml	1.8		
Ø x L in mm	-		
Cat. No.	cryo tubes		
rotor Cat. No. 2426-B			
Cat. No.	-		
boring Ø x L in mm	12.5 x 36		
tubes per rotor	20		
max. RCF	21,130		
radius in mm	84		
run-up in sec	19		
run-down in sec, braked	28		
temperature in °C1)	+4		

capacity in ml	0.2	0.2		
Ø x L in mm	6 x 18	-		
Cat. No.	-	PCR strips		
	<b>Q</b>			
rotor Cat. No. <b>2418-A</b>				
Cat. No.	-			
boring Ø x L in mm	6.5 x 15.5			
tubes per rotor	32 4 x 8			
max. RCF	14,338			
radius in mm	57			
run-up in sec	19			
run-down in sec, braked	28			
temperature in °C1)	+4			

#### Optional for rotor 2424-B, 2426-B, 2427-A and 2430-B



lid **Cat. No. 2425**with bio-containment<sup>3)</sup>,
autoclavable



lid **Cat. No. 2423**with bio-containment <sup>3)</sup>,
autoclavable,
phenol-resistant

#### Optional for rotor 2418-A



lid Cat. No. E3243

 $<sup>^{3)}\,</sup>$  Tested by the TÜV in conformity with DIN EN 61010, section 2-020.

TECHNOLOGY	MIKR	0 200	MIKRO 200 R		
Microlitre centrifuge, without rotor	clas	ssic	cooled		
Power supply*)	200 – 240 V 1 ~ 100 – 127 V 1 ~		200 – 240 V 1 ~	100 – 127 V 1 ~	
Frequency	50 – 6	60 Hz	50 Hz	50 – 60 Hz	
Consumption	240 VA	270 VA	400 VA	630 VA	
Emission, Immunity	EN / IEC 61326-1, class B	FCC class B	EN / IEC 61326-1, class B	FCC class B	
Max. capacity	30 x 1.5 / 2.0 ml				
Max. RPM (speed)	15,000 min <sup>-1</sup>				
Max. RCF	21,382				
Running time	1 sec - 99 min : 59 sec, ∞ continuous run, short cycle mode (impulse key)				
Dimensions (HxWxD)	260 x 275	x 344 mm	260 x 281 x 547 mm		
Weight	approx.	11.5 kg	approx. 28 kg		
Refrigeration					
Temperature control, infinitely variable	-		from -10 to +40 °C		
Cat. No.	2400	2400-01	2405	2405-01	

<sup>\*)</sup> Other voltages on request.







Hettich centrifuges comply with all relevant EU standards in effect and conform to the European level of quality and safety for medical devices. Evidence is provided by national and international test marks such as IEC 61010 or the CE conformity. The ISO 9001, ISO 13485 and ISO 14001 certificates accredited to the company bear witness to the extreme care and responsibility Hettich puts into the manufacturing of centrifuges and their accessories.



Föhrenstr. 12 D-78532 Tuttlingen Germany www.hettichlab.com

www.hettichlab.com info@hettichlab.com service@hettichlab.com Phone +49 (0)7461/705 -0 Fax +49 (0)7461/705 -122

National Sales: -200
International Sales: -201
National Service: -202
International Service: -203

